



EIE-1-TH

Endur ID Infant Easy Thermal Identification

Robert Chadwick
Endur ID Inc.
8 Merrill Industrial Dr., Unit 4
Hampton, NH 03842
PH: (603) 758-1488 EXT 102
FX: (603) 758-1118

bchadwick@endurid.com
www.endurid.com



NICU & INFANTS

Identification Challenges

- ▶ **Small Size – difficulty accommodating information**
- ▶ **Skin Sensitivity – not “soft” enough**
- ▶ **Security**
- ▶ **Cleaning and Durability**
- ▶ **Difficult to accommodate all situations and sizes**
- ▶ **The need for tethers to reach outside the swaddle**
- ▶ **Not easy to use**
- ▶ **Secure closure – adhesive issues**



Endur ID EIE-1-TH Infant / NICU Band





Endur ID EIE-1-TH Infant / NICU Band

- ▶ Very Soft Silicone Material to prevent skin irritation
- ▶ Clasp Closure for security and tamper evident



- ▶ Smooth on the “Baby” side for comfort
- ▶ Printed with Standard Thermal Printer - Sharp text and barcodes
- ▶ Large Print Area
- ▶ Can be cleaned with Alcohol
- ▶ Can be used 6 different ways and long tethers are a possible option
- ▶ Will fit all size infants and small children
- ▶ No PVC, Latex or Phalates



Option 1 Tail On Top






Option 2 Infant Easy, Pass Through Hole






Option 3

Infant Easy, Pass Through Hole Double Back Tail






Option 4 Infant Tether





Option 5 NICU





Option 6 NICU Tether

